



## HHW CARE PRODUCTS (I) PVT LTD

Our firms are formerly known as  
1. Handicapped Care Products India Pvt. Ltd., 2. Dinesh Industrial Corporation, 3. Lotus Cycle and Rickshaw Co.

### FOLDING WHEELCHAIR IN CHROME FINISH WITH MAG WHEEL MAZIC

8" CASTOR WITH P.U.TUBELESS TYRE

24' GLASS-FILLED NYLON WHEEL WITH PLASTIC RUSH RIM  
AND P.U. TUBLESS TYRES

FOOTREST-P.P MOLDED WITH ZINC PLATED HEIGHT  
ADJUSTABLE LEG REST

FRAME- M.S.SILVER CHROME COATED

PLASTIC BUCKET WITH LID

SEAT & BACKSEAT- REMOVABLE & WASHABLE

ARMREST-FIXED WITH SOFT P.U.TOP

BOTH SIDE HAND OPERATED ZINC PLATED BRAKES

SEALED BALL BEARINGS ARE USED IN WHEEL FOR A  
SMOOTH RIDE

**Weight 15.0kg**  
**Seat width 16" & 18"**  
**Load Capacity 100kg**

**Price - 7030/-**



**Office :**  
584, Bhagirathpura , Indore (M.P.) 452003  
Ph.: 0731-2423404,2422260

**Factory :**  
79, SDA Compound, Dewas Naka,  
Indore (M.P.) 452010. Ph.: 0731 4023053

Email : [info@activeforall.co.in](mailto:info@activeforall.co.in) ♦ Website : [www.activeforall.co.in](http://www.activeforall.co.in)